

## PRODUCT DESCRIPTION

Id-Fungi Plates is a culture medium with a transparent membrane, allowing growth of moulds, yeasts & dermatophytes.

### PRINCIPLE:

Antibiotic and pH are specifically settled to promote fungi growth. The membrane is limiting moulds growth inside the medium but ensure a good transfer of nutrients.

### PRESENTATION-APPEARANCE:

Medium ready to use: box of 20 plates

### STORAGE:

At reception, store Id-fungi plates between 4 & 8°C.

### LAPSING:

If storage conditions are respected, Id-Fungi Plates are useable for 3 months after production.

### IN WHICH CASE USE THE ID-FP:

#### a. Purification :

With an inoculation loop, sample gently the colony of interest and settled it to the center of the plate. You can also carve a plug with a scalpel at the younger limit of the colony and settled it to the center of the plate.

#### b. Culture of the sample (first intention) :

Dispose the sample directly to the center of the plate.

*Alternative 1:* Make an exhaustion streak of your sample carefully on the membrane.

*Alternative 2:* Dispose your sample at different spots onto the membrane.

#### c. Transplanting :

With a sterile dry swab or an inoculation loop take the sample – you can also use a scalpel to carve a plug - and put it to the center of the plate.

Incubation between 25°C & 37 °C directly linked to the sample nature.

Reading results between 24h & 96h linked to the identification method used (MALDI-TOF/ Micro-Macroscopic observations) and to the species.

## PERFORMANCES & QUALITY CONTROL OF ID-FUNGI PLATES

### Fertility evaluation:

- *Candida albicans* ATCC 10231 48-72h 20-25°C
- *Aspergillus brasiliensis* ATCC 16404 48-72h 20-25°C
- *Trichophyton rubrum* IP 1464.83 5J 37°C

### Bacterial growth inhibition:

- *Escherichia coli* ATCC 8739 48-72h 25°C

### Id-FP performances were tested with two kinds of samples.

#### 1. Collection Fungi species.

More than 21 clinical or environmental species have been tested:

*Penicillium* (2), *Aspergillus* (4), *Fusarium* (1), *Alternaria* (2), *Cladosporium*(2), *Trichoderma* (1), *Mucorals* (2), *Stachybotris* (1), yeasts (2), Dermatophytes (5), ...

**Environmental strains:** Growth visible from 24h incubation at 25°C.

**Yeasts:** growth visible from 24h incubation 25°C.

**Dermatophytes:** growth visible from 48-72h incubation at 25°C.

#### 2. Patient samples.

In progress

### Wastes treatment

Remove all materials (biologic or non-biologic) as the local regulation required.

### BIBLIOGRAPHIC REFERENCES

In progress